

HUGS TRAINING TEST A

Name: _____ Date: _____

- T or F** When a patient leaves the hospital, it is necessary to discharge the Hugs tag from the system (at the PC) before it is removed from the patient.
- T or F** The Hugs tag must be removed before giving a patient a bath.
- T or F** The rounded part of the Hugs tag should face the patient's foot.
- T or F** The Hugs tag should always be removed before the patient leaves the safe area.
- T or F** A "heartbeat" is sent by the Hugs tag to ensure that it is working properly.
- T or F** The Hugs tag can be soaked for cleaning.
- T or F** To save time, it is OK to put the strap in both ears of the Hugs tag (make a loop) and later slip it on the patient's ankle.
- T or F** The Hugs System should be explained to parents.
- T or F** If a Hugs tag is too tight, and a strap needs to be changed, you must discharge the tag and get both a new tag and a new strap.

How tight should a Hugs tag be? _____

How often should you assess the Hugs tag? _____

When is the "Transport" function used? _____

When will a "Tamper Alarm" sound? _____

What does the "BabySense™" feature do? _____
