



PEDZ SOFTWARE USAGE QUIZ

Objective:

The purpose of these exercises is to affirm your familiarization with the Pedz system software and to check how you utilize it. Carefully read each question; and give a brief, clear answer in the space provided.

1. You put a Pedz tag on a patient, the Smith girl in room 314. **A:** *Where do you enter the patient information on the computer?* **B:** *Clearly write out the patient information as you would type it into the computer:*
 - a) _____
 - b) _____

2. Patient Smith is moving from room 314 to room 360. **A:** *What needs to be done at the computer?*
 - a) _____

3. Patient Smith needs to leave the floor to get some tests. You will be leaving the tag on her. She will be out of the department for 45 minutes. **A:** *What needs to be done at the computer?*
 - a) _____

4. Patient Smith has returned from her tests. **A:** *What needs to be done at the computer?*
 - a) _____

5. Child Smith's parents think the Pedz band is too tight. You cannot pull the band back because that will set off a tamper alarm; you cannot just cut the band either. **A:** *What needs to be done at the computer before you remove the old band and put a new band on?*
 - a) _____

6. Child Smith is ready to go home. **A:** *What needs to be done at the computer before you remove the tag?*
 - a) _____

7. A TAMPER alarm is underway, and it has been determined that the child is secure. The child's father thought the tag/band was on too tight and pulled back on the tag! Before you can clear the alarm, the computer asks for a note about the alarm. **A:** *Clearly write out the appropriate note you might type into the computer.*
 - a) _____